

<b>WIRED/WIRELESS</b>	
<b>PROGRAMS CONFIGURATION</b>	
<b>CHANNELS</b>	Number of channels of stimulation
<b>PROGRAM CATEGORIES</b>	<b>TYPE OF PROGRAMS</b>
<b>REHABILITATION I</b>	Atrophy, Reinforcement, Motor point...
<b>PAIN RELIEF I</b>	Tens Gate Control, Endorphins...
<b>VASCULAR</b>	Heavy legs, Capillarization...
<b>CONDITIONING I</b>	Resistance, Strength, Recovery...
<b>REHABILITATION II</b>	ACL, Rotator Cuff, Lumbar Stabilization...
<b>PAIN RELIEF II</b>	Low back Pain, Cervical Pain, Epicondylitis...
<b>NEURO REHABILITATION</b>	Spacticity, Hemiplegic Foot & Shoulder...
<b>CONDITIONING II</b>	Endurance, Explosive Strength, Massage...
<b>INCONTINENCE</b>	Stress, Urge, Mixed, Post Partum Prevention...
<b>DENERVATED</b>	Partial and Total, Automatic and Manual...
<b>DIRECT CURRENT</b>	Iontophoresis, Edema, Hyperhidrosis...
<b>ADDITIONAL FUNCTIONS</b>	
<b>2+2 OR +TENS OPTION</b>	Combine 2 different programs
<b>MOTOR POINT PEN</b>	Localize muscle motor point for optimal electrode placement
<b>SKIP</b>	Move directly to the next phase
<b>PROGRAMMABILITY</b>	Create your own programs
<b>MI-SCAN</b>	Customize the stimulation for comfort and efficiency
<b>MI-TENS</b>	Work efficiently behind the muscle threshold
<b>MI-RANGE</b>	Work efficiently at good intensity for low frequency programs
<b>MI-ACTION</b>	Active sessions with voluntary contraction detected by the device
<b>MULTI-SESSIONS™</b>	2 remotes + Smart docking station to manage several sessions independently
<b>TRIGGER</b>	Starts voluntary working phase
<b>SOFTWARE</b>	Create your patient and history lists
<b>AUTOMATIC RAMPE SEARCH</b>	Partial denervation
<b>WARRANTY</b>	3 years



REHAB	THETA	PHYSIO	WIRELESS PRO 2CH STANDARD	WIRELESS PRO 4CH STANDARD	WIRELESS PRO 4CH FULL	
Wired	Wired	Wired	Wireless	Wireless	Wireless	
Standard	Full	Optimum	Standard	Standard	Full	
4	4	4	2 Modules	4 Modules	4 Modules	
X	X	X	X	X	X	STANDARD
X	X	X	X	X	X	
X	X	X	X	X	X	
X	X	X	X	X	X	
	X	X			X	FULL
	X	X			X	
	X	X			X	
		X				OPTIMUM
		X				
		X				
X	X	X	X	X	X	
X	X	X	X	X	X	
X	X	X	X	X	X	
X	X	X	X	X	X	
Mi-ready	2 Mi-cables	4 Mi-cables	2 Ch	4 Ch	4 Ch	
Mi-ready	2 Mi-cables	4 Mi-cables	2 Ch	4 Ch	4 Ch	
Mi-ready	2 Mi-cables	4 Mi-cables	2 Ch	4 Ch	4 Ch	
Mi-ready	2 Mi-cables	4 Mi-cables	2 Ch	4 Ch	4 Ch	
				X	X	
	X	X			X	
			X	X	X	
		X				

**WIRELESS PRO** is today's most advanced electrotherapy technology offering enhanced ease-of-use and convenience while saving space and time to optimize patient treatment. It offers a new way of stimulation without constraints allowing you and your patient a total freedom of movement facilitating patient manipulation and combined exercises.

#### **MULTI-SESSIONS Concept**

Discover the MULTI-SESSIONS technology offered by WIRELESS PRO four channel devices. Four modules, two remotes and a smart docking station allow you to synchronize modules and remotes to independently manage several sessions simultaneously. Additional remotes/modules are optional.